

Human Survivin Promoter (6309 bp) Luciferase Reporter Construct

CH3 Catalog # 110100 **Unit Size:** 1 µg

PRODUCT OVERVIEW	Human Survivin promoter (5' Sal I – Hind IIIA) inserted into luciferase reporter construct based upon a plasmid backbone similar to pGEM3.
PURIFICATION	Standard mini-preparation from Amp ^r culture.
RECOMMENDED USAGE	See plasmid map for restriction enzyme sites. Transfection of mammalian cells using the method of choice of the laboratory investigator.
STORAGE	Store plasmid DNA at -20° to -80° or bacterial cells at -80° until expansion of transformed bacteria for scaled-up plasmid preparation of vector DNA prior to transfection studies.

All products are for Research Use Only and not for in vitro diagnostic use


This Package Insert is provided for product evaluation purposes only and is not intended to be used in place of the Package Insert shipped with the product.

University Commons ♦ Suite 8 ♦ 1416 Sweet Home Road ♦ Amherst, NY 14228 ♦ 716.688.5222

www.ch3biosystems.com

Human Survivin Promoter (6309 bp) Luciferase Reporter Construct

CH3 Catalog # 110100 Unit Size: 1 µg



```

SalI Hind III
1  GTCGACAAGC TTCCTCTGCA CCCGAGAGTG GTGTGCATGT CCCAAGAAGC
51  CAGCTATATC CATCCACCAC CCACCCCACC ACACCGCCTC TTTCTCTCAC
101 ACCGAATAGC ATTCAAGAAT TTCCACTGTG GGCTGGGCGC AGGGGCTAGT
151 GCCTGTAATC CCAGCACTTT AGGAGGCCGA GGTGGGCGAA TCACCTGAGG
201 TCAGGAGTTC AAGACCAGCC TGGCCAACAT GGCAAAACCC TGTCTCTACT
251 AAAAAATACAA AAATTAGCCA GGCGTGGTGG CTGGCGCCTG TAATCCTAGA
301 TACGTGGGAG GCTGAGTTAG GAGAATCTCT TGAATCCGGG AGGCAGAGGT
351 TGCAGTGAGC TGAGATTGTG CCCCTGCACT CCAGCCTGGG TGACAAAGCC
... ..
6001 AAGTGTGGG ATTACAGGCG TGAGCCACTG CACCCGGCCT GCACGCGTTC
6051 TTTGAAAGCA GTCGAGGGGG CGCTAGGTGT GGCAGGGAC GAGCTGGCGC
6101 GCGTCGCTG GGTGCACCCG GACCACGGGC AGAGCCACGC GCGGGGAGGA
6151 CTACAACCTC CGGCACACCC CGCGCCGCC CGCCTCTACT CCCAGAAGGC
6201 CGCGGGGGT GGACCGCCTA AGAGGGCGTG CGCTCCCGAC ATGCCCCGGG
6251 GCGCGCCATT AACC GCCAGA TTTGAATCGC GGGACCCGTT GGCAGAGGTG
6301 GCGCGGCGG CAAAGCTT
Hind III
    
```

All products are for Research Use Only and not for in vitro diagnostic use

This Package Insert is provided for product evaluation purposes only and is not intended to be used in place of the Package Insert shipped with the product.

University Commons ♦ Suite 8 ♦ 1416 Sweet Home Road ♦ Amherst, NY 14228 ♦ 716.688.5222

www.ch3biosystems.com