


SURVIVIN Mammalian Expression Vector pcDNA3-HA-Survivin T34A

CH3 Catalog # 131005 Unit Size: 1 µg



	<i>HindIII</i>	HA Tag		<i>Human Survivin</i>
1	AAGCTT	TATGT	ATCCGTATGA TGTTCCTGAT TATGCTGGTG	CCCCGACGTT
51	GCCCCCTGCC	TGGCAGCCCT	TTCTCAAGGA CCACCGCATC	TCTACATCA
				Mutation (A->G)
101	AGAACTGGCC	CTTCTTGAG	GGCTGCGCCT GCGCCCCGGA	GCGGATGGCC
151	GAGGCTGGCT	TCATCCACTG	CCCCACTGAG AACGAGCCAG	ACTTGGCCCA
201	GTGTTTCTTC	TGCTTCAAGG	AGCTGGAAGG CTGGGAGCCA	GATGACGACC
251	CCATAGAGGA	ACATAAAAAG	CATTCGTCCG GTTGCCTTT	CCTTCTGTC
301	AAGAAGCAGT	TTGAAGAATT	AACCCTTGGT GAATTTTGA	AACTGGACAG
351	AGAAAGAGCC	AAGAACAAAA	TTGCAAAGGA AACCAACAAT	AAGAAGAAAG
401	AATTTGAGGA	AACTGCGAAG	AAAGTGCGCC GTGCCATCGA	GCAGCTGGCT
451	GCCATGGATT	GAGGCCTCTG	GCCGGAGCTG CCTGGTCCCA	GAGTGGCTGC
		BamHI		
501	ACCACTTCCG	GATCC		

All products are for Research Use Only and not for in vitro diagnostic use

This Package Insert is provided for product evaluation purposes only and is not intended to be used in place of the Package Insert shipped with the product.

University Commons ♦ Suite 3 ♦ 1416 Sweet Home Road ♦ Amherst, NY 14228 ♦ 716.688.5222

www.ch3biosystems.com