

## **SURVIVIN Bacterial Expression Vectors For Survivin Research**

**CH3 Catalog Series # 132000 Unit Size: 1 µg**

<b>PRODUCT OVERVIEW</b>	Human survivin inserted into pGEX4T3 for production of GST fusion proteins.
<b>PURIFICATION</b>	Standard mini-preparation from Amp <sup>r</sup> culture.
<b>RECOMMENDED USAGE</b>	Transfection of bacterial cells using the method of choice of the laboratory investigator.
<b>STORAGE</b>	Store plasmid DNA at -20° to -80° or bacterial cells at -80° until expansion of transformed bacteria for scaled-up plasmid preparation of vector DNA prior to transfection studies.

*All products are for Research Use Only and not for in vitro diagnostic use*

*This Package Insert is provided for product evaluation purposes only and is not intended to be used in place of the Package Insert shipped with the product.*

---

University Commons ♦ Suite 3 ♦ 1416 Sweet Home Road ♦ Amherst, NY 14228 ♦ 716.688.5222

[www.ch3biosystems.com](http://www.ch3biosystems.com)