

SURVIVIN shRNA Targeting Vectors For Survivin Research

CH3 Catalog Series # 140000 Unit Size: 1 µg

PRODUCT OVERVIEW	BIRC5 (survivin) inserted into retroviral vector V2HS_94582/3 with or without EGFP, luciferase).
PURIFICATION	Standard mini-preparation from Amp ^r culture.
RECOMMENDED USAGE	Transfection of mammalian cells using the method of choice of the laboratory investigator.
STORAGE	Store plasmid DNA at -20° to -80° or bacterial cells at -80° until expansion of transformed bacteria for scaled-up plasmid preparation of vector DNA prior to transfection studies.

All products are for Research Use Only and not for in vitro diagnostic use

This Package Insert is provided for product evaluation purposes only and is not intended to be used in place of the Package Insert shipped with the product.

University Commons ♦ Suite 3 ♦ 1416 Sweet Home Road ♦ Amherst, NY 14228 ♦ 716.688.5222

www.ch3biosystems.com