

## Survivin ΔEx3-Specific Antibody

CH3 Catalog # 150002 Unit Size: 100 µl

<b>PRODUCT OVERVIEW</b>	Anti-Survivin ΔEx3 is an IgY-purified chicken polyclonal antibody that was raised against the survivin Δ3-specific peptide (RKPTIRRKNLRKLRK).
<b>PURIFICATION</b>	Affinity purification by agarose chromatography using conjugated immunizing peptide.
<b>FORMULATION</b>	PBS/glycerol/0.02% sodium azide.
<b>RECOMMENDED USAGE</b>	Western Blotting.
<b>TESTING</b>	Tested using protein extracts from HeLa cells transfected either with survivin, survivin 2B or survivin ΔEx3. This reagent recognizes only survivin ΔEx3. For experimental procedures, consult Li, F. and Ling, X. 2006 ( <a href="http://www.ncbi.nlm.nih.gov/pubmed/16557517">http://www.ncbi.nlm.nih.gov/pubmed/16557517</a> ).
<b>STORAGE</b>	Avoid freeze/thaw cycles. Store at 4° or prepare appropriately sized aliquots and keep at -20° C.

*All products are for Research Use Only and not for in vitro diagnostic use*

*This Package Insert is provided for product evaluation purposes only and is not intended to be used in place of the Package Insert shipped with the product.*

---

University Commons ♦ Suite 8 ♦ 1416 Sweet Home Road ♦ Amherst, NY 14228 ♦ 716.688.5222

[www.ch3biosystems.com](http://www.ch3biosystems.com)