

Survivin 3B-Specific Antibody

CH3 Catalog # 150004 Unit Size: 100 µl

PRODUCT OVERVIEW	Anti-Survivin 3B is an IgY-purified chicken polyclonal antibody that was raised against survivin 3B peptide (CZKNKIERALLAE).
PURIFICATION	Affinity purification by agarose chromatography using conjugated immunizing peptide.
FORMULATION	PBS/glycerol/0.02% sodium azide.
RECOMMENDED USAGE	Western Blotting.
TESTING	Tested using protein extracts from HeLa cells transfected with either survivin, survivin 2B or survivin ΔEx3. This reagent recognizes only survivin 3B. For experimental procedures, consult Li, F. and Ling, X. 2006 (http://www.ncbi.nlm.nih.gov/pubmed/16557517).
STORAGE	Avoid freeze/thaw cycles. Store at 4° or prepare appropriately sized aliquots and keep at -20° C.

All products are for Research Use Only and not for in vitro diagnostic use

This Package Insert is provided for product evaluation purposes only and is not intended to be used in place of the Package Insert shipped with the product.

University Commons ♦ Suite 8 ♦ 1416 Sweet Home Road ♦ Amherst, NY 14228 ♦ 716.688.5222

www.ch3biosystems.com