

HepG2 Cell Line

CH3 Catalog # 930001 **Unit Size:** Flask or Vial

PRODUCT OVERVIEW	Human hepatocellular carcinoma cell line.
PURIFICATION	NA
FORMULATION	25 to 50 X 10 ⁶ cells in culture flask containing RPMI 1640 culture medium and 5% fetal bovine serum or frozen cell pellet in cryo-protectant.
RECOMMENDED USAGE	Cell culture. Useful to provide cell proteins for large-scale purification in downstream applications. Cells adapt well to a wide variety of culture media. Fetal bovine serum at 5% is recommended.
TESTING	Cultured on tissue culture plastic in 5% fetal bovine serum with either DMEM or RPMI 1640 culture media.
STORAGE	Place in culture immediately upon arrival; 35° humidified CO ₂ incubator. Use 5 to 7.5% CO ₂ .

For Research Use Only. Not for in vitro diagnostic use.