

LUCIFERASE REPORTER CONSTRUCTS For Survivin Reserach

CH3 Catalog Series # 110000 Human #120000 Mouse

Unit Size: 1 µg

PRODUCT OVERVIEW

Luciferase reporter inserted into pGEM3-like plasmid backbone. For general map see #110000.

#110000 Promoter-less

#110001 Promoter-less with Neo cassette

#11000x With human survivin promoter

#12000x With mouse survivin promoter

PURIFICATION

Standard mini-preparation from Amp^r culture.

RECOMMENDED USAGE

Transfection of mammalian cells using the method of choice of the laboratory investigator.

STORAGE

Store plasmid DNA at -20° to -80° or bacterial cells at -80° until expansion of transformed bacteria for scaled-up plasmid preparation of vector DNA prior to transfection studies.

All products are for Research Use Only and not for in vitro diagnostic use

This Package Insert is provided for product evaluation purposes only and is not intended to be used in place of the Package Insert shipped with the product.

University Commons ♦ Suite 3 ♦ 1416 Sweet Home Road ♦ Amherst, NY 14228 ♦ 716.688.5222

www.ch3biosystems.com