

## **LNCaP Cell Line**

**CH3 Catalog # 940002 Unit Size:** Flask or Vial

<b>PRODUCT OVERVIEW</b>	Human prostate adenocarcinoma cell line.
<b>PURIFICATION</b>	NA
<b>FORMULATION</b>	10 to 20 X 10 <sup>6</sup> cells in flask containing RPMI 1640 culture medium and 5% fetal bovine serum or frozen cell pellet in cryo-protectant. Maintained as described by <a href="#">Horoszewicz JS et al., Cancer Res. 43:1809-18, 1983.</a>
<b>RECOMMENDED USAGE</b>	Cell culture. Useful for large-scale cell culture to provide human cellular proteins for purification and analysis in downstream applications. Adapt well to a wide variety of culture media. Fetal bovine serum at 5% is recommended.
<b>TESTING</b>	Cultured on tissue culture plastic in 5% fetal bovine serum with RPMI 1640.
<b>STORAGE</b>	Place in culture immediately upon arrival; 35-37° humidified CO <sub>2</sub> incubator. Use 5 to 7.5% CO <sub>2</sub> .

*All products are for Research Use Only and not for in vitro diagnostic use  
This Package Insert is provided for product evaluation purposes only and is not intended to be used in place of the Package Insert shipped with the product.*