



MARCH 2019

**PREMIUM PACKAGE --- CH3 BIOSYSTEMS CATALOG #940001:
E0771 MURINE BREAST CANCER CELL LINE**

Selected Lots of the E0771 cell line are authenticated by short tandem repeat profiling and screened for the following pathogens:

(M. spp., M. pulmonis)	Mycoplasma
(TMEV)	Theiler's Murine Encephalomyelitis Virus
(PVM)	Pneumonia Virus of Mice
(MVM)	Minute Virus of Mice
(MHV)	Mouse Hepatitis Virus
(MPV1-5)	Mouse Parvovirus
(C. bovis)	Corynebacterium bovis
	Sendai Virus

In addition, species-specific PCR evaluation is performed to rule out any contamination by rat, human, hamster or monkey cells.

*Screening for additional pathogens may be available upon request.
Contact us for customized pricing.*