

HeLa Cell Line

CH3 Catalog # 940003 Unit Size: Flask/Vial

PRODUCT OVERVIEW	Human cervical cancer cell line.
PURIFICATION	NA
FORMULATION	10 to 20 x 10 ⁶ cells in flask containing RPMI 1640 culture medium and 5% fetal bovine serum or frozen cell pellet in cryoprotectant.
RECOMMENDED USAGE	Cell culture. Useful for large-scale cell culture to provide human cellular proteins for purification and analysis in downstream applications. Adapt well to a wide variety of culture media. Fetal bovine serum at 5% is recommended.
TESTING	Cultured on tissue culture plastic in 5% fetal bovine serum with either DMEM or RPMI 1640.
STORAGE	Place in culture immediately upon arrival; 35-37° humidified CO ₂ incubator. Use 5% CO ₂ .

*All products are for Research Use Only and not for in vitro diagnostic use
This Package Insert is provided for product evaluation purposes only.*

University Commons ♦ Suite 3 ♦ 1416 Sweet Home Road ♦ Amherst, NY 14228 ♦ 716.688.5222

www.ch3biosystems.com