

## LNCaP Cell Line

**CH3 Catalog # 940002 Unit Size:** Flask/Vial

<b>PRODUCT OVERVIEW</b>	Human prostate adenocarcinoma cell line.
<b>PURIFICATION</b>	NA
<b>FORMULATION</b>	10 to 20 x 10 <sup>6</sup> cells in flask containing RPMI 1640 culture medium and 5% fetal bovine serum or frozen cell pellet in cryo-protectant. Maintained as described by <a href="#">Horoszewicz JS et al., Cancer Res. 43:1809-18, 1983.</a>
<b>RECOMMENDED USAGE</b>	Cell culture. Useful for large-scale cell culture to provide human cellular proteins for purification and analysis in downstream applications. Adapt well to a wide variety of culture media. Fetal bovine serum at 5% is recommended.
<b>TESTING</b>	Cultured on tissue culture plastic in 5% fetal bovine serum with RPMI 1640.
<b>STORAGE</b>	Place in culture immediately upon arrival; 35-37° humidified CO <sub>2</sub> incubator. Use 5% CO <sub>2</sub> .

*All products are for Research Use Only and not for in vitro diagnostic use  
This Package Insert is provided for product evaluation purposes only.*